

Universal Library™

Redistribution Guide



**MEASUREMENT
COMPUTING™**

Document Revision 2.8, June, 2008

© Copyright 2008, Measurement Computing Corporation

Your new Measurement Computing product comes with a fantastic extra —

Management committed to your satisfaction!

Refer to www.mccdaq.com/execteam.html for the names, titles, and contact information of each key executive at Measurement Computing.

Thank you for choosing a Measurement Computing product—and congratulations! You own the finest, and you can now enjoy the protection of the most comprehensive warranties and unmatched phone tech support. It's the embodiment of our mission:

- To provide PC-based data acquisition hardware and software that will save time and save money.

Simple installations minimize the time between setting up your system and actually making measurements. We offer quick and simple access to outstanding live FREE technical support to help integrate MCC products into a DAQ system.

Limited Lifetime Warranty: Most MCC products are covered by a limited lifetime warranty against defects in materials or workmanship for the life of the product, to the original purchaser, unless otherwise noted. Any products found to be defective in material or workmanship will be repaired, replaced with same or similar device, or refunded at MCC's discretion. For specific information, please refer to the terms and conditions of sale.

Harsh Environment Warranty® Program: Any Measurement Computing product that is damaged due to misuse, or any reason, may be eligible for replacement with the same or similar device for 50% of the current list price. I/O boards face some harsh environments, some harsher than the boards are designed to withstand. Contact MCC to determine your product's eligibility for this program

30 Day Money-Back Guarantee: Any Measurement Computing Corporation product may be returned within 30 days of purchase for a full refund of the price paid for the product being returned. If you are not satisfied, or chose the wrong product by mistake, you do not have to keep it.

These warranties are in lieu of all other warranties, expressed or implied, including any implied warranty of merchantability or fitness for a particular application. The remedies provided herein are the buyer's sole and exclusive remedies. Neither Measurement Computing Corporation, nor its employees shall be liable for any direct or indirect, special, incidental or consequential damage arising from the use of its products, even if Measurement Computing Corporation has been notified in advance of the possibility of such damages.

Licensing Information

Each original copy of Universal Library is licensed for development use on one CPU at a time. It is theft to make copies of this program for simultaneous program development. If a customer creates an application using the Universal Library, they may distribute the necessary runtime files (Universal Library driver files) with their application royalty free. They may not distribute any files that give their customer the ability to develop applications using the Universal Library.

Trademark and Copyright Information

TracerDAQ, Universal Library, Harsh Environment Warranty, Measurement Computing Corporation, and the Measurement Computing logo are either trademarks or registered trademarks of Measurement Computing Corporation.

Windows, Microsoft, and Visual Studio are either trademarks or registered trademarks of Microsoft Corporation

LabVIEW is a trademark of National Instruments.

CompactFlash is a registered trademark of SanDisk Corporation.

XBee and XBee-PRO are trademarks of MaxStream, Inc.

All other trademarks are the property of their respective owners.

Information furnished by Measurement Computing Corporation is believed to be accurate and reliable. However, no responsibility is assumed by Measurement Computing Corporation neither for its use; nor for any infringements of patents or other rights of third parties, which may result from its use. No license is granted by implication or otherwise under any patent or copyrights of Measurement Computing Corporation.

All rights reserved. No part of this publication may be reproduced, stored in a retrieval system, or transmitted, in any form by any means, electronic, mechanical, by photocopying, recording, or otherwise without the prior written permission of Measurement Computing Corporation.

Notice

Measurement Computing Corporation does not authorize any Measurement Computing Corporation product for use in life support systems and/or devices without prior written consent from Measurement Computing Corporation. Life support devices/systems are devices or systems which, a) are intended for surgical implantation into the body, or b) support or sustain life and whose failure to perform can be reasonably expected to result in injury. Measurement Computing Corporation products are not designed with the components required, and are not subject to the testing required to ensure a level of reliability suitable for the treatment and diagnosis of people.

Table of Contents

Overview	5
Installing the Universal Library for Windows Vista, Windows XP, and Windows 2000.....	5
Installing <i>InstaCal</i> with the Universal Library	7
Running a Silent <i>InstaCal</i> Installation	9
Installing Universal Library support for the .NET Framework.....	10
Tips and techniques.....	10
Digital signing.....	10
Deploying applications with devices that have serial numbers.....	10
Deploying applications when using PCI boards	10

Redistributing Applications with the Universal Library

Overview

The information in this document explains how to distribute applications developed using the 32-bit Universal Library (UL). Developers wishing to redistribute the files associated with the Universal Library and *InstaCal* may do so only if they are not altered or edited in any way.

This document includes information on how to install the Measurement Computing device drivers, Universal Library, and *InstaCal* under Windows Vista, Windows XP, and Windows 2000.

These instructions assume that Measurement Computing-specific files are installed to a directory named C:\Program Files\Measurement Computing\DAQ.

Installing the Universal Library for Windows Vista, Windows XP, and Windows 2000

This section describes explains how to install the Measurement Computing device drivers, the Universal Library, and *InstaCal* under Windows Vista, Windows XP, and Windows 2000. The general steps to install an application are:

1. Copy your application files to a directory on your system.
2. Update the registry and system files as needed by your application.
3. Create a desktop folder for your program and associated files.
4. Restart the computer, if necessary.

These instructions describe how to install the Universal Library manually, and to supply an appropriate CB.CFG for the installed Measurement Computing hardware. These instructions assume that *InstaCal* has not been installed. For instructions on installing *InstaCal*, refer to ""Installing *InstaCal* with the Universal Library" on page 7.

In addition to your application files, the following Universal Library files must be installed:

File description	File name
Universal Library Windows 2000, Windows XP, and Windows Vista device drivers	<ul style="list-style-type: none">▪ CBUL32.sys▪ CBULWDM.sys▪ pcidaqlib.sys▪ mccusb.sys▪ usbdqlib.sys▪ usblddaqlib.sys
Universal Library DLL files	<ul style="list-style-type: none">▪ CBW32.DLL▪ DaqLib.DLL▪ mcsktsifc.dll
Measurement Computing UL support files	<ul style="list-style-type: none">▪ CBERRCODE.TXT▪ ULPROPS.TXT
Hardware configuration files	<ul style="list-style-type: none">▪ CB.CFG (created using <i>InstaCal</i>)▪ CBI95.INF▪ CBICOM.INF (For Measurement Computing COM port devices only)▪ DAQLIB.INF▪ MCCUSB.INF

In addition, the Windows system configuration must be updated to access the Universal Library.


To update the system configuration:

1. Copy CBW32.DLL, CBERCODE.TXT, ULPROPS.TXT, and DaqLib.DLL to a directory on your hard drive, such as C:\Program Files\Measurement Computing\DAQ.

2. Update the Windows registry to include the following key:

[HKEY_LOCAL_MACHINE]\Software\Universal Library\RootDir

3. Add a New String Value to the key that was just created: RootDir = C:\Program Files\Measurement Computing\DAQ.



4. Copy CB.CFG (created using *InstaCal*) to a directory on your hard drive:

- When running Windows Vista, copy CB.CFG to C:\ProgramData\Measurement Computing\DAQ.
- When running on Windows platforms other than Windows Vista, copy CB.CFG to C:\Program Files\Measurement Computing\DAQ.


Note: Steps 5 and 6 apply to Windows Vista only. If you are running Windows XP or Windows 2000, proceed to step 7.

5. When running Windows Vista, update the Windows registry to include the following key:

[HKEY_LOCAL_MACHINE]\Software\Universal Library\ConfigDir


6. When running Windows Vista, add a New String Value to the new ConfigDir key:

ConfigDir = C:\ProgramData\Measurement Computing\DAQ. This value sets the directory where the CB.CFG file is created when you run *InstaCal*.



7. Update the PATH environment variable to include the directory containing CBW32.DLL: Edit the following registry key as shown below:

[HKEY_LOCAL_MACHINE]\System\CurrentControlSet\Control\Session Manager\Environment



8. Install the device driver files: copy CBUL32.sys, CBULWDM.sys, pcidaqlib.sys, mccusb.sys, usbdqlib.sys, and usblddaqlib.sys to the System Drivers directory (usually C:\WINDOWS\system32\drivers).
9. Update the Windows registry to include the following keys:


```
[HKEY_LOCAL_MACHINE]\System\CurrentControlSet\Services\CBUL32\CBUL32
```

```
[HKEY_LOCAL_MACHINE]\System\CurrentControlSet\Services\CBUL32
```

Add the following fields to the [HKEY_LOCAL_MACHINE]\System\CurrentControlSet \Services\CBUL32 key to include the device driver (cbul32.sys). Refer to the image below:

```
DisplayName = Measurement Computing DataAcq
ImagePath = System32\Drivers\CBUL32.SYS
Type = 1
Start = 1
Group = Extended Base
ErrorControl = 1
```

Note that Type, Start and ErrorControl must all be DWORD values, while the rest are String values.



10. Copy CBI95.INF, CBICOM.INF, and DAQLIB.INF, and MCCUSB.INF to the Windows INF directory (usually C:\WINDOWS\inf).

Note

If you install a Measurement Computing plug-and-play device before you install the device driver, you must find the device in the device manager and update the driver. Refer to the ReadMe.txt file for details.

Installing InstaCal with the Universal Library

The following instructions are applicable to Windows Vista, Windows XP, and Windows 2000.

The CB.CFG configuration file must be available to the Universal Library, and must contain information that is accurate with respect to the installed Measurement Computing devices. We recommend that you run *InstaCal* on each target machine to update the CB.CFG when the hardware configuration has changed.

InstaCal also provides field-calibration and testing for Measurement Computing devices. *InstaCal* requires the Universal Library DLL and device driver to be installed, as explained in "Installing the Universal Library for Windows Vista, Windows XP, and Windows 2000" on page 5. To distribute *InstaCal* with your application, add the following files to the distribution:


File description	File name
<i>InstaCal</i> executable	INSCAL32.EXE

File description	File name
InstaCal support files	<ul style="list-style-type: none"> LOGO.BMP NDERCODE.TXT SCANVIEW.EXE CBI_CAL.DLL CBI_NODE.DLL CBI_PROP.DLL CBI_TEST.DLL VICOMPONENTS.OCX
System support files	<ul style="list-style-type: none"> MFC42.DLL MSVCRT.DLL OLEPRO32.DLL OLEAUT32.DLL <p>It is likely that these files already exist on the target machine. Check the version of existing files against those being installed. More recent versions on the target machine should not be overwritten. These files should be installed to the system directory – usually C:\Windows\System32</p>

InstaCal requires the following changes to the Windows registry:

1. Copy the *InstaCal* files *inscal32.exe*, *logo.bmp*, *ndercode.txt*, *scanview.exe*, *cbi_cal.dll*, *cbi_node.dll*, *cbi_prop.dll*, *cbi_test.dll*, and *vicomponents.ocx* to the C:\Program Files\Measurement Computing\DAQ directory.
2. Create the registry key:

[HKEY_LOCAL_MACHINE]\Software\ComputerBoards



3. Copy the system support files *mfc42.dll*, *olepro32.dll* and *oleaut32.dll* to the system directory – usually C:\Windows\System32.

Note

The support files *mfc42.dll*, *olepro32.dll*, and *oleaut32.dll* may already exist on the target machine. Do not overwrite them if they are a more recent version than those you are trying to install.

4. The support files *MFC42.DLL*, *OLEPRO32.DLL*, *OLEAUT32.DLL*, and *VICOMPONENTS.OCX* must also be registered. If they already exist on the target machine, they may be in use (and therefore locked). The full installation and registration for some or all of these files may only occur at boot time.

The files can be made to self-register when Windows restarts. Create the following keys under:

[HKEY_LOCAL_MACHINE]\Software\Microsoft\Windows\CurrentVersion\RunOnce


with these values:

CBI1 = C:\Windows\System32\regsvr32 /s C:\Windows\System32\mfc42.dll

CBI2 = C:\Windows\System32\regsvr32 /s C:\Windows\System32\olepro32.dll

CBI3 = C:\Windows\System32\regsvr32 /s C:\Windows\System32\oleaut32.dll

CBI4 = C:\Windows\System32\regsvr32 /s C:\Program Files\Measurement Computing\DAQ\vicomponents.ocx



Distribution software also provides solutions to handle the complexity of this situation, providing the capacity to:

- Overwrite only older versions of these files.
- Delay overwriting until restarting Windows.
- Allow the files to self-register upon restart.
- Prompt the user to restart the computer.

The necessary support files will register on startup.

Running a Silent *InstaCal* Installation

Perform a "Silent Install" to use the *InstaCal* installer instead of copying the files and editing the registry yourself. Your installer will launch the *InstaCal* installation, but you will not see the installation window or any prompts from the installation when using the `"/qn "` parameter.

To include the *InstaCal* setup program, do the following:

1. Add the following to your installation code to launch the file:

```
msiexec /i "path to InstaCal.msi" /qn REBOOT=ReallySuppress
```

2. Place all *InstaCal* files in a folder:

▪ CBI95.IN	▪ logo.bmp
▪ CBICOM.IN	▪ MCCBoardCfg.mde
▪ DAQLIB.IN	▪ ReadMe.txt
▪ ICalUL.msi	▪ ulprops.txt

When performing a silent installation, *InstaCal* installs to the Measurement Computing\DAQ directory in the Program Files directory (usually C:\Program Files\Measurement Computing\DAQ).

The installation will not prompt the user to reboot the computer. Since rebooting is required, the master installation should either prompt the user to reboot or force them to do so.

The CB.CFG configuration file is created the first time that you run *InstaCal*.

- When running Windows Vista, the CB.CFG file is created in C:\ProgramData\Measurement Computing\DAQ. You must be logged into *InstaCal* with Administrative rights in order to make configuration changes.
- When running Windows XP and Windows 2000, the CB.CFG file is created in C:\Program Files\Measurement Computing\DAQ.

Installing Universal Library support for the .NET Framework

The .NET support is an extension of the Universal Library. It requires that both the Microsoft .NET Framework and *InstaCal* or Universal Library be installed. Refer to the previous sections for information about installing *InstaCal* and/or the Universal Library.

To install Universal Library support for the Microsoft .NET Framework, do the following:

1. Install either *InstaCal* or the Universal Library.
2. Copy the .NET Universal Library file *MccDaq.DLL* to your target installation directory.

When *MccDaq.DLL* is in the same directory as your application, you do not need to add it to the Global Assembly Cache.

Tips and techniques

This section includes tips and suggestions regarding common problems that are encountered when deploying an application. This list is intended only to provide guidance to those who are redistributing their applications with the Universal Library, and may not apply to all situations.

Digital signing

In order to avoid "Unknown Publisher" messages when installing a custom application, you should digitally sign the installation driver and application package with your unique Software Certificate. This will notify your customers that you are the author of the software being installed.


Deploying applications with devices that have serial numbers

If the application you are deploying is to be used with a device that has a custom serial number — such as a USB or Ethernet based device — you may need to alter that configuration file to allow it to operate with the device regardless of the serial number.

Note

You only need to alter the configuration file when you install the Universal Library with your installation and are providing your own *CB.CFG* configuration file. If you are installing *InstaCal* and intend on the user running *InstaCal*, you do not need to alter the configuration file.

When a device has a serial number, the configuration file stores the serial number as miscellaneous option 16 (Misc Option [16]). If your device does not have a serial number that matches this configuration item, such as when you run *InstaCal* or replace the device), an error will occur. To avoid this error and force the Universal Library to ignore the serial number of the device, change miscellaneous option 16 to "0" in the configuration file:



Deploying applications when using PCI boards

Devices that are assigned resources dynamically—such as PCI boards—require the use of *InstaCal* for board configuration to determine the system resources assigned to the board.

Measurement Computing Corporation
10 Commerce Way
Suite 1008
Norton, Massachusetts 02766
(508) 946-5100
Fax: (508) 946-9500
E-mail: info@mccdaq.com
www.mccdaq.com